


Zero BINDING™ Staining Destaining Trays for Gels & Blots


Special Surface Re-usable Trays with Lids Designed for Electrophoresis & Blotting

Product Highlights

- Multi-purpose trays ideal for mini gels and blots
- Requires low volume of staining & destaining solutions
- Lid protects evaporation and inhalation of solvents
- Minimum surface absorption of dye
- Minimum absorption of protein and antibodies
- Applicable for almost all Western blots
- Applicable for Northern and Southern blot assays
- Deep design for faster shaking of gels and blots
- Corners designed for easy decantation of liquids
- Re-usable for several times
- Microwave safe



Applications




Application for gel staining and destaining on shakers or rocker without evaporation



Useful for gel bands observation and storage in liquid medium



Viewing and photography of gel under gel light box illuminator



NRL Zero Binding Trays are useful for storage of gels over years in water. Two mini gels (14% SDS-TG) stained with NRL silver stain was stored in Zero Binding Trays for more than 2 years.

Ordering Information

Catalog #	Product Description	Price, \$
SDT-001	Zero BINDING™ Staining Trays (pack of 5)	80.00
SDT-002	Zero BINDING™ Staining Trays (pack of 10)	150.00